
Instructions for use	Breathing circuit set, dual limb	4
Gebrauchsanweisung	Beatmungsschlauchset, mit zwei Schenkeln.....	7
Instructions d'utilisation	Kit de circuit respiratoire à deux branches	10
Instrucciones de uso	Equipo respiratorio de dos ramas	13
Istruzioni per l'uso	Set circuito paziente, branca doppia	16
Instruções de uso	Kit de respiração, alça dupla	19
使用说明	呼吸装置, 双肢.....	22
使用説明書	デュアルリム呼吸回路セット.....	24
Инструкции по эксплуатации	Комплект дыхательного контура с двумя патрубками.....	27
Kullanım talimatları	Solunum devresi seti, çift parçalı.....	30


REF 260241, 260244

CE0197 **ERC**

HAMILTON
MEDICAL

1 To patient
Zum Patienten
Port inspiratoire
Hacia el paciente
Al paziente
Para o paciente
至病人
吸気側接続口
Порт «к пациенту»
Hastaya

4 Inspiratory limb
Inspirationsschenkel
Branche inspiratoire
Rama inspiratoria
Branca inspiratoria
Alça inspiratória
吸気肢
吸気回路
Патрубок вдоха
Inspiratuvar parça

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Conventions

⚠ WARNING

Alerts the user to the possibility of injury, death, or other serious adverse reactions associated with the use or misuse of the device.

⚠ CAUTION

Alerts the user to the possibility of a problem with the device associated with its use or misuse, such as device malfunction, device failure, damage to the device, or damage to other property.

NOTICE

Emphasizes information of particular importance.

1 Overview

The dual limb breathing circuit set (PN 260241, 260244) can be used with all Hamilton Medical ventilators.

1.1 Intended use

The neonatal/pediatric breathing circuit sets, with PN 260241 and PN 260244, comprise breathing circuits that include a Y-piece, connectors, and elbow adapters. Each of these breathing sets is intended to be used to connect the HAMILTON-C1/C2/C3/C6/T1/MR1/G5/S1 during invasive or non invasive ventilation of neonatal or pediatric patients.

The breathing set is a medical device intended for use by qualified, trained personnel.

2 Safety information

⚠ WARNING

- Reprocessing of Hamilton Medical single-use products can affect the product properties and may cause injury to the patient. For example, a change to the surface structure during reprocessing may lead to a change in the tear strength or cause actual cracking. Furthermore, an altered surface structure may result in, for example, a microbial aggregation of spores, allergens, and

pyrogens, or cause an increase in the number of particles released as a result of chemical changes in the material properties.

- Hamilton Medical does not assume any warranty for the proper function of single-use items if they are reprocessed and reused by the user. Any liability on the part of Hamilton Medical for any damages resulting from the use of reprocessed or reused single-use parts, or use of any non- Hamilton Medical accessories and consumables, is excluded.
- Follow hospital infection control procedures for all tests and use.
- Remove the safety cover(s) before use.
- Discard the product if it fails two subsequent preoperational checks or if there is any sign of damage.
- To prevent patient or ventilator contamination, always use a bacteria filter or HMEF between the patient and the inspiratory port. If no bacteria filter is used, the exhaled gas can contaminate the ventilator.

- Adding attachments or other parts/ assemblies to a breathing system can change the pressure gradient across the ventilator, which can adversely affect ventilator performance.
- Connect only a single breathing circuit at a time. Do not connect any circuit as an extension to another one. Doing so can seriously harm or kill the patient.
- Route the breathing circuit without tension and without any kinks from the ventilator or humidifier to the patient, and protect against unintended influences.

⚠ CAUTION

- *(USA only) Federal law restricts this device to sale by or on the order of a physician.*
- *Handle used products as contaminated goods according to local laws and regulations or hospital internal procedures.*
- *Always carry out all preoperational checks before connecting the patient to prevent possible patient injury. Make sure another source of ventilator support is available.*

NOTICE

- Before use, read the appropriate ventilator *Operator's Manual*.
- Replace the product in accordance with hospital infection control procedures, or depending on the patient's secretions and nebulization of medication. Consider using a filter at the expiratory port during nebulization.
- Any incident with the device leading to serious patient injury, death, or a potential threat to public health must be reported to the manufacturer and the relevant authorities.
- Any bacteria filter, HMEF, or additional accessories in the expiratory limb may substantially increase flow resistance and impair ventilation.
- When adding parts to the Hamilton Medical breathing circuit configurations, do not exceed the inspiratory and expiratory resistance values of the ventilator breathing system as specified in the ventilator *Operator's Manual*.
- Certain parts of the breathing circuit may be visible on the MRI image. Keep the breathing circuit out of the scanned area as much as possible.

3 Connecting and calibrating components

To connect the breathing circuit

- ▶ Connect the breathing circuit as shown on page 2, and ensure that the limbs are connected to the appropriate ports on the ventilator.

4 Calibration and preoperational check

Perform the preoperational check before patient use, including calibrating the flow sensor or pressure line and performing a tightness test. Follow the instructions provided in the ventilator *Operator's Manual*.

5 Delivery

The components are delivered clean and ready for clinical use.

6 Disposal

A used component must be handled as contaminated. Follow all local, state, and federal regulations with respect to waste management and environmental protection when disposing of used parts.

7 Environmental requirements and storage conditions

Environmental requirements and storage conditions

Operating temperature -15°C to 40°C (5°F to 104°F)
5% to 95% relative humidity, noncondensing

Storage temperature -20°C to 60°C (-4°F to 140°F)
5% to 95% relative humidity, noncondensing

8 Technical specifications

Technical data	PN 260241 (1.50 m)	PN 260244 (3.00 m)
----------------	--------------------	--------------------

Compliance (at 60 cmH ₂ O)	0.3 ml/cmH ₂ O	0.6 ml/cmH ₂ O
---------------------------------------	---------------------------	---------------------------

Conforms to	EN 12342 / ISO 5367	
-------------	---------------------	--

Tube length	1.50 m	3.00 m
-------------	--------	--------

Conical tube connectors	Conforms to ISO 5356-1	
-------------------------	------------------------	--

Smallest internal diameter	8.0 mm	
----------------------------	--------	--

Max. operating duration	28 days	
-------------------------	---------	--

Maximum working pressure	200 cmH ₂ O	
--------------------------	------------------------	--

Resistance to flow at	2.5 l/min	15 l/min
-----------------------	-----------	----------

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Specifications are subject to change without notice. For all proprietary trademarks (®) and third-party trademarks (§) used by Hamilton Medical AG, see www.hamilton-medical.com/trademarks. © 2020 Hamilton Medical AG. All rights reserved.
Hamilton Medical AG
medin Medical Innovations GmbH