Instructions for use
Gebrauchsanweisung
Instrucciones de uso

IntelliCuff

REF 951001

Software version 1.0.x.x
624741/06 | 2017-05-18

CE0197

**HAMILTON
MEDICAL**
Intelligent Ventilation since 1983

IntelliCuff

Instructions for use

April 2017

624741/06

1	English	5
2	Deutsch	21
3	Español	39

Table of Contents

1	Device overview	8
2	Getting started	9
	Insert the batteries and connect the USB cable	9
	Connect the device to primary power	9
	Turn on the device	9
	Connect the cuff pressure tube to the device and patient	9
	Check for airway and cuff leaks	10
	Review and adjust pressure setting	10
	Lock/unlock the control panel (optional)	10
	Attach the device to the mount	10
3	Performing a time-limited hold	11
4	Deflation and extubation	11
	4.1 Deflating the cuff	12
	4.2 Turning off the device	12
5	Powering the device	12
6	Alarms and troubleshooting	13
7	Configuration	14
	7.1 Configuring units of measurement and default target pressure	14
	7.2 Configuring time-limited hold settings	14
8	Mounting the IntelliCuff controller	15
	8.1 Attaching IntelliCuff to a fixed mount	15
	8.2 Attaching IntelliCuff to the hanger	16
9	Cleaning, maintenance, and disposal	16
	9.1 Cleaning	16
	9.2 Maintenance	16
	9.3 Disposal	16
10	Intended use and audience	16
	10.1 Intended use	16
	10.2 Intended audience	17
11	Standards and approvals	17
12	Specifications	17
	12.1 Physical, performance, and environmental data	17
	12.2 Symbols on the device labels	18
	12.3 System overview	18
	12.4 Essential performance	19
13	Parts and accessories	19
14	Warranty	19
	14.1 Miscellaneous	19

© 2017 Hamilton Medical AG. All rights reserved. Printed in Switzerland.

No part of this publication may be reproduced, stored in a database or retrieval system, or transmitted in any form or by any means, electronic, mechanical, or by photocopying, recording, or otherwise, without prior written permission of Hamilton Medical AG.

This document may be revised, replaced, or made obsolete by other documents by Hamilton Medical AG at any time and without notice. Ensure that you have the most current applicable version of this document; if in doubt, contact the technical support department of Hamilton Medical AG, Switzerland. While the information set forth herein is believed to be accurate, it is not a substitute for the exercise of professional judgment.

Nothing in this document shall limit or restrict in any way Hamilton Medical AG's right to revise or otherwise change or modify the equipment (including its software) described herein, without notice. In the absence of an express, written agreement to the contrary, Hamilton Medical AG has no obligation to furnish any such revisions, changes, or modifications to the owner or user of the equipment (including software) described herein.

The equipment must be operated, serviced, or upgraded only by trained professionals. Hamilton Medical AG's sole responsibility with respect to the equipment and its use is as stated in the limited warranty provided in the operator's manual.

Hamilton Medical AG shall not be liable for any loss, cost, expense, inconvenience, or damage that may arise out of misuse of the product, or if non-Hamilton Medical AG parts were used when replacing parts, or if serial numbers were amended, deleted, or removed.

If returning parts to Hamilton Medical AG, be sure to use the standard Hamilton returned goods authorization (RGA) procedure. Disposal of parts shall follow all local, state, and federal regulation with respect to environmental protection.

Product and company names mentioned herein may be the trademarks and/or registered trademarks of their respective owners.

Hamilton Medical AG will make available, on request, circuit diagrams, component parts lists, descriptions, calibration instructions, or other information that will assist appropriately trained personnel to repair those parts of the equipment designated by Hamilton Medical AG to be repairable.

Document conventions

WARNING

A WARNING alerts the user to the possibility of injury, death, or other serious adverse reactions associated with the use or misuse of the device.

CAUTION

A CAUTION alerts the user to the possibility of a problem with the device associated with its use or misuse, such as device malfunction, device failure, damage to the device, or damage to other property.

NOTICE

A NOTE emphasizes information of particular importance.

Units of measure

The document uses cmH₂O representative of all pressure units. 1 cmH₂O equals 0.981 mbar, which equals 0.981 hPa.

General notes

WARNING

- Additional equipment connected to medical electrical equipment must comply with the respective IEC or ISO standards (for example, IEC 60950 for data processing equipment). Furthermore, all configurations shall comply with the requirements for medical electrical systems (see IEC 60601-1, clause 16). Anybody connecting additional equipment to medical electrical equipment is configuring a medical system and is, therefore, responsible that the system complies with the requirements for medical systems.
- MR UNSAFE. Keep away from magnetic resonance imaging (MRI) equipment. IntelliCuff poses unacceptable risks to patient, medical staff, or other persons within the MR environment.
- Modifications to the device are not permitted.
- To prevent increased emissions, decreased immunity, or interrupted operation of the IntelliCuff device or any accessories, use only accessories or cables that are expressly stated in this manual.

CAUTION

- *Use only Hamilton Medical disposable tubing with a filter and safety valve. Use of any other tubing will result in the immediate loss of cuff pressure if disconnected on the ventilator end. Use of any other tubing without a filter may result in the device being contaminated.*
- *DO NOT kink the tubing.*

NOTICE


- The use of this equipment is restricted to one patient at a time who is intubated with an endotracheal tube or tracheal cannula (both are referred to as *ET tube* in this document).
- If there is visible damage to any part of the IntelliCuff device, do not use the device. Technical service is required.
- Familiarize yourself with these *Instructions for use* before using this device on a patient.
- Do not simultaneously touch conductive components (for example, the USB port) or conductive parts of the cuff pressure controller, the enclosure, and the patient.
- Install the device in a position where the primary power supply can easily be disconnected.

- To electrically isolate the IntelliCuff from all poles of the primary power supply simultaneously, disconnect the power plug.
- The device is not protected against the effects of defibrillator use.
- The manufacturer can only be responsible for the safety, reliability, and performance of the IntelliCuff device if all of the following requirements are met:
 - Appropriately trained personnel carry out assembly operations, extensions, readjustments, modifications, maintenance, or repairs.
 - The electrical installation of the relevant room complies with the appropriate requirements.
 - The IntelliCuff device is used in accordance with the *Instructions for use*.
- Portable and mobile RF communications equipment can affect the IntelliCuff device.
- The IntelliCuff device requires special precautions regarding EMC, and must be installed and put into service according to the EMC information provided in the *EMC declarations* (PN 624750).

1 Device overview


Not all elements shown here are displayed at the same time. They are shown for illustration purposes only.


Alarm and power display

- 1 Alarm silence button
- 2 Alarm lamp
- 3 Alarm icons
- 4 Power On/Off button
- 5 Power icons


Pressure-related displays and controls

- 1 Pressure trend (increase/decrease) indicator
- 2 Icons related to time-limited hold
- 3 Displayed when target pressure can be adjusted
- 4 Current pressure
- 5 Unit of measure
- 6 Target (set) pressure
- 7 Adjustment buttons (+/-)

Control panel

- 1 Decrease (-) target pressure
- 2 Increase (+) target pressure
- 3 Deflate cuff
- 4 Lock/unlock control panel
- 5 Perform time-limited hold


IntelliCuff connection

- 1 USB power cable
- 2 IntelliCuff device
- 3 Cuff pressure tube with filter
- 4 Cuff connector with ET tube


2 Getting started

Note that the device accepts your changes after 3 seconds of inaction, and confirms the change with a beep.

Insert the batteries and connect the USB cable

1. Lift up the cover using a screwdriver.
2. Insert freshly charged batteries, as shown on the battery case.
3. Connect the mini-USB connector end of the power cable to the port.
4. Close the cover.

Connect the device to primary power

- ▶ Connect the USB power cable to any of the following:
 - USB connector on the ventilator
 - Power adapter; then plug the adapter into an AC power source

The display shows the power indicator and the battery charge level. See Section 5.

WARNING

Do not remove the batteries without the device being connected to primary power supply. The device will cease operation.

Turn on the device

- ▶ Long-press the **Power On/Off** button.

The device performs a self test, during which all of the display lights are lit and the buzzer sounds. If this does not occur, see Section 6.

When the target pressure on the display blinks, the self test is successful and the device starts applying the set default pressure. After 3 seconds of inaction, the blinking stops.

Note that all alarms are silenced for 1 minute after the self test to provide time to connect the patient.


Connect the cuff pressure tube to the device and patient




- ▶ If not already connected, connect the cuff pressure tube as shown. See page 8.

Once connected, the device starts applying the set target pressure.

NOTICE

- Only use a Hamilton Medical cuff pressure tube with filter (PN 282016).
- ET tube with high-volume cuffs may take longer to inflate.

Check for airway and cuff leaks

- **Leaks in the airway.** IntelliCuff may help to reduce airway leakages around the cuff. You may increase the cuff pressure based on leakages detected by the ventilator. For pressure settings above 30 cmH₂O, consider using a larger ET tube. Also check the patient's throat for bubbling or gurgling sounds.
- **Leaks in the cuff.** When present, the cuff-leakage alarm is generated on the IntelliCuff device.

Review and adjust pressure setting

For an overview of the pressure display, see Section 1.

As shown here, the device is increasing the pressure, currently at 5 cmH₂O, to the target of 25 cmH₂O.

NOTICE

We recommend that you keep cuff pressure (default 25 cmH₂O) below 30 cmH₂O for adult/pediatric ET tubes, and below 20 cmH₂O for neonatal ET tubes.

To adjust the target pressure

1. Press the Increase (+) or Decrease (-) button to enable adjustment of the pressure.
The target pressure blinks and the word SET appears.
2. Press the Increase (+) or Decrease (-) button while the value is blinking and set the pressure as desired.

After 3 seconds of inaction, the numbers stop blinking, the device beeps, and the new setting is applied.

Lock/unlock the control panel (optional)

- ▶ Long-press the **Lock** button to lock or unlock the keypad.
When locked, the Lock icon is displayed. All buttons are disabled except for the **Alarm silence** button.

Attach the device to the mount

See Section 8 for more details.


3 Performing a time-limited hold

NOTICE

Any changes made to the hold pressure or time are reset to the factory defaults once the device is turned off.

The maximum pressure allowed under any circumstances is limited to a total of 55 cmH₂O.

Figure 1. Time-limited hold controls and icons

- 1 Hold button
- 2 The set target plus the hold pressure
- 3 Hold countdown timer (1 dot = 1 minute)
- 4 Hold icon

The hold function temporarily increases the cuff pressure by a specified amount for a set period of time to more completely seal the airway and prevent aspiration.

By default, the hold is active for 5 minutes, and applies 5 cmH₂O above the currently set pressure.

You can change the increased pressure setting in 5 cmH₂O increments from a minimum of 5 cmH₂O to a maximum of 25 cmH₂O, and set the hold duration to either 5 or 10 minutes.

For details on changing the duration of the hold, see Section 7.

To perform a time-limited hold

1. Long-press the **Hold** button for > 3 seconds. The device beeps when Hold mode is initiated, and the Hold icon appears.

The target pressure in the bottom half of the display, showing the target pressure plus the configured addition, blinks for 3 seconds.

2. Press the Increase (+) or Decrease (-) button while the value is blinking and set the pressure as desired.

To accept the pressure setting, do nothing. After 3 seconds, the displayed value stops blinking and the Hold maneuver begins.

During the Hold maneuver

The pressure increases to the set amount.

The hold timer, comprising a Watch icon and countdown dots, is shown on the display. Each dot represents one minute. When set to 5 minutes, 5 dots plus the Watch icon are displayed; when set to 10 minutes, 10 dots are displayed.

As the time counts down, the watch icon and the current minute dot flash continuously. As each minute passes, the dots disappear one by one, from right to left until the hold is complete.

At the end of the configured time

- The device beeps
- The pressure returns to the pre-hold target pressure
- The Hold icon and countdown timer disappear

To stop the hold at any time

- ▶ Touch the **Hold** button again on the control panel to cancel the hold and return operation to the previously set target pressure.

4 Deflation and extubation

After deflating the cuff, you can extubate the patient and turn off the device.

Figure 2 provides an overview of the deflation controls and icons.

Figure 2. Cuff deflation controls and icons

- 1 Deflate button
- 2 Target pressure set to 0
- 3 Pressure trend arrow (decreasing)
- 4 Deflate icon
- 5 Cuff deflated for > 1 minute alarm

4.1 Deflating the cuff

To deflate the cuff

1. Long-press the **Deflate** button for > 3 seconds.
The device beeps and displays the Deflate icon to indicate that deflation is in progress.
The target pressure display shows a blinking 0 for 3 seconds. Then the device beeps and applies negative pressure until the current pressure is 0. The decrease pressure trend arrow is displayed.
When the pressure reaches 0, the current and target pressure both display 0. The device beeps again to indicate that the cuff is deflated.
2. Extubate the patient.
3. Disconnect the device from the ET tube.
4. Turn off the device within 1 minute of deflating the cuff.
If the cuff is deflated for > 1 minute without being either reinflated (*Deactivation impossible* alarm) or the device being turned off, you are reminded that you deflated the cuff and that the cuff is still in the deflated state. See Table 2.

4.2 Turning off the device

NOTICE

You cannot turn off the device when the pressure measurement is > 5 cmH₂O. In this case, a *Deactivation impossible* alarm is generated. See Table 2.

To turn off the device

1. Deflate the cuff fully and ensure both the current and target pressure values are 0.
2. Long-press the **Power On/Off** button for > 3 s.

5 Powering the device

NOTICE

- Battery life indications are approximate. The actual battery life depends on settings, battery age, and level of battery charge. To ensure maximum battery life, maintain a full charge and minimize the number of complete discharges.
- Replace batteries once their capacity drops below acceptable levels.

The device is powered by two rechargeable AA batteries and a USB cable connected to AC power.

Two new, fully charged batteries provide approximately 5 hours of standard operation.

The device display indicates the connected power supplies, as described in Table 1. See also page 8.

Table 1. Power indicators

	Power supply icon, visible when a primary power supply is connected
	67% to 100% charge level, approximately 3 to 5 hours operating time
	34% to 66% charge level, approximately 1 to 3 hours operating time
	10% to 33% charge level, approximately 1 hour operating time
	0% to 9% charge level, no more than 0.5 hours operation; connect to a primary power supply immediately

6 Alarms and troubleshooting

When an alarm is generated, the device emits audible beeps and the alarm lamp lights yellow or red, depending on the alarm priority. The IntelliCuff device has two alarm priorities: high and medium. See Tables 2 and 3 for details.

To silence an alarm

- ▶ Review the alarm, and if appropriate, press the **Alarm silence** button.
The alarm is silenced for 2 minutes and the Alarm silence icon appears (Section 1).

Table 2. IntelliCuff alarms

Alarm type	Alarm lamp	Audio response	Action required
High-priority alarm	Red, flashing	A sequence of beeps, repeated until the alarm is reset.	Depends on the alarm; see Table 3
Medium-priority alarm	Yellow, flashing	A sequence of beeps, repeated periodically.	Depends on the alarm; see Table 3
Technical fault ^{1,2}	Red, flashing	A sequence of beeps, repeated until the device is turned off, if technically possible. At a minimum, a continuous buzzer tone; the buzzer cannot be silenced.	<ul style="list-style-type: none"> • Exchange the device • Contact Hamilton Medical for a replacement
Loss of power ^{1,2}	Alarm lamp not functional without power	Continuous buzzer tone for 1 minute; the buzzer cannot be silenced.	<ul style="list-style-type: none"> • Connect the device to primary power supply • Replace the batteries

1. All LCD segments blink. The number in the lower display area denotes the code of the technical failure.
2. Device is in safe state. Cuff remains inflated, but is not regulated anymore.

Table 3. Troubleshooting alarms

Alarm name	Alarm type	Possible causes	Action required
 Cuff system leakage	High-priority alarm	<ul style="list-style-type: none"> • Cuff loses pressure • The cuff pressure tube is not correctly connected 	<ul style="list-style-type: none"> • Check applied pressure • Check cuff pressure tube, ET tubing, all connections • Change ET tube, if needed
 Battery low	High-priority alarm	The batteries are depleted	<ul style="list-style-type: none"> • Make sure the device is connected to a primary power source • Replace the batteries
 Battery 10%	Medium-priority alarm	Ten percent or less of battery power remains	<ul style="list-style-type: none"> • Connect the device to a primary power source • Replace the batteries
Pressure above set limit	Medium-priority alarm	The device is unable to meet the specified pressure setting	<ul style="list-style-type: none"> • Check cuff pressure tube, ET tubing, all connections • Replace the cuff tubing
 Cuff deflated	Medium-priority alarm	The cuff has been deflated for over 1 minute	Disconnect the cuff and turn off the device, or re-inflate the cuff
 Deactivation impossible Current pressure blinks	Medium-priority alarm	Cuff pressure is above 5 cmH ₂ O when trying to turn off the device	Deflate or disconnect the ET-tube cuff before turning off the device

Alarm name	Alarm type	Possible causes	Action required
 Technical fault All display lights blink	Technical fault	Various	<ul style="list-style-type: none"> Remove the device from use Contact Hamilton Medical for a replacement
Loss of power ¹	Loss of power	Primary power supply disconnected or failed, and the batteries are depleted	<ul style="list-style-type: none"> Make sure the device is connected to a primary power source Replace the batteries

1. The device resumes operation with the previous settings once a power supply has been reestablished.

Should multiple alarms occur at the same time, all corresponding alarm icons are lit and the alarm lamp flashes in the color of the highest-priority alarm.

5. Once set, do nothing for 3 seconds to exit Configuration.
The device beeps and the display stops blinking. Normal operation continues.

7 Configuration

You can set units of measurement and default target pressure, and change the hold settings.

7.1 Configuring units of measurement and default target pressure

To access general Configuration mode

The device must be turned off when you start.

1. Simultaneously long-press the **Power On/Off** and **Alarm silence** buttons for > 3 seconds to start the device in Configuration mode.

The controller beeps, performs a self-test, and displays blinking units of measure to show it has entered Configuration mode.

2. Use the Increase (+), and Decrease (-) buttons to move left or right to select the desired unit: hPa, mbar, or cmH2O.

3. Once set, do nothing for 3 seconds.

The unit stops blinking and cycles to pressure settings, displaying a blinking target pressure.

4. Use the Increase (+) and Decrease (-) buttons to change the default target pressure.

NOTICE

If you accidentally hold down the Increase (+) button, the pressure setting increases to a maximum of 50. The device beeps to alert you as long as you hold down the button.

If you hold down the Decrease (-) button, the pressure setting decreases to a minimum of 5. The device beeps to alert you as long as you hold down the button.

7.2 Configuring time-limited hold settings

Any changes you make to the time-limited hold settings apply only until the device is turned off. When it is turned on again, the device resets to the factory default values.

To change the time-limited hold settings

1. Long-press the **Hold** button > 7 seconds.

The device beeps, the Watch icon appears, the + icon next to the target pressure appears, and the target pressure blinks.

2. Use the Increase (+) and Decrease (-) buttons to change the amount above the target pressure to apply during the hold. Set in increments of 5, from a minimum of 5 to a maximum of 25.

3. After setting the value, either do nothing for 3 seconds or press the **Hold** button again to cycle to the next setting, **Hold duration**.

The display value changes to minutes, and the **Min** unit appears.

4. Use the Increase (+) and Decrease (-) buttons to change the time. Options are 5 or 10.

5. Do nothing for 3 seconds.

The device automatically exits the Configuration mode and starts the hold maneuver, increasing the pressure to the new target value.

8 Mounting the IntelliCuff controller

Several mounting options are available:

- Fixed mounts (Multi mount, Rail mount)
- A hanger for IntelliCuff

The back of the IntelliCuff controller is designed to mount directly to the attachment clip provided with both mount types.

8.1 Attaching IntelliCuff to a fixed mount

A fixed mount allows attachment of the controller to a rail.

Figure 3. IntelliCuff attachment well (1)



Figure 4. Fixed mount (Multi mount)

- 1 Attachment clip
- 2 Latch to attach/release device
- 3 Tab

To attach the controller to a fixed mount


1. Holding the controller with one hand, fit the top of the attachment well (1 in Figure 3) on the back of the controller over the tab at the top of the attachment clip (3 in Figure 4).

The tab fits under the upper lip of the attachment well (1 in Figure 5).

2. Using the index or middle finger of the hand holding the controller, push the attach/release

latch forward (toward you), and at the same time, push the device onto the clip (2 and 3 below).

Figure 5. Attach to mount


The controller connects to the clip with an audible click.

3. Pull gently on the controller to ensure it is firmly attached.

To detach the device from the mount


- ▶ Grasping the attached controller with one hand, push the attachment/release latch toward the controller, while at the same time pulling the controller away and up from the clip, over the tab (3 in Figure 4).

8.2 Attaching IntelliCuff to the hanger

The hanger allows you to hang the controller as needed.

Figure 6. IntelliCuff hanger


- 1 Attachment clip
- 2 Latch to attach/release device

To attach the controller to the hanger

1. Holding the controller in one hand, and the hanger in the other, fit the top of the attachment well (1 in Figure 3) on the back of the controller over the tab at the top of the attachment clip on the hanger (1 in Figure 6).
2. Push the attach/release latch forward, and at the same time, push the controller and the clip together.
The controller connects to the clip with an audible click.
3. Pull gently on the controller to ensure it is firmly attached.

To detach the device from the hanger

1. Holding the hanger in one hand and the attached controller in the other face down, push the attachment/release latch toward the controller, while at the same time pushing the upper portion of the handle against the controller.
This helps lever the clip on the hanger up and away from the controller.

2. Once the clip is released, slide the hanger down and out of the attachment well.

9 Cleaning, maintenance, and disposal

9.1 Cleaning

⚠ WARNING

Always disconnect the device from electrical power before cleaning.

NOTICE

- Strong solvents, such as acetone or trichloroethylene, may damage the surface.
- Be sure to only clean around the connection ports, not inside them.
- Be particularly careful with infectious patients, and follow your hospital infection protocol procedures.

Use a lint-free cloth dampened with any of the following:

- Warm water (maximum 40°C (104°F)) and soap
- A dilute and nonacid agent
- A surfactant
- A cleaning agent in a base of ammonia or alcohol

9.2 Maintenance

Preventive maintenance is not mandatory, but if preventive maintenance is required by hospital protocol, Hamilton Medical recommends to perform the tests once per year as described in the service documentation.

9.3 Disposal

Dispose of all parts removed from the device according to your institution's protocol. Follow all local, state, and federal regulations with respect to environmental protection, especially when disposing of the electronic device or parts of it, such as batteries.

10 Intended use and audience

10.1 Intended use

The IntelliCuff device is intended to continuously measure and automatically maintain the user-set cuff pressure of an endotracheal tube (ETT) or tracheostomy tube (TT) during mechanical ventilation.

The device can be used with any mechanical ventilator, as follows:

- When used with a non-Hamilton Medical ventilator, IntelliCuff adjusts the cuff pressure to values set on the device.
- When used with a Hamilton-C6 ventilator, IntelliCuff adjusts the cuff pressure to values set either on the device or on the ventilator, depending on configuration.

The device is to be used during ventilation of adults, pediatrics, and neonates, who are intubated with ETT or TT, in the following areas:

- In the intensive care ward or in the recovery room
- In the operation room during intubation narcosis
- For emergency medical care or primary care
- During transport within and outside of the hospital
- During transfer by rescue vehicles, ship, jet, or helicopter

CAUTION

(USA only): Federal law restricts this device to sale by or on the order of a physician.

10.2 Intended audience

This IntelliCuff device is a medical device intended for use by qualified, trained personnel under the direction of a licensed physician and within the limits of its stated technical specifications.

Read and understand these *Instructions for use* to complete your training on the use of this device.

11 Standards and approvals

IntelliCuff was developed in accordance with pertinent international standards and FDA guidelines.

The device is manufactured within the standards specified by EN ISO 13485 and EN ISO 9001, Council Directive 93/42/EEC, Annex II, Article 1 certified quality management system.

The device meets the essential requirements of Council Directive 93/42/EEC. It is a class IIa device.

The device meets relevant parts of, among others, the following standards:

IEC 60601-1:2012 Medical electrical equipment
IEC 60601-1-2:2007 EMC
IEC 60601-1-6:2013 + IEC 62366:2014 Usability
IEC 60601-1-8:2012 Alarm standard
IEC 60601-1-12:2014 Medical emergency equipment
EN 13718-1:2008 Medical devices in air ambulances
EN 1789:2010 Medical devices in road ambulance

12 Specifications

12.1 Physical, performance, and environmental data

Table 4. Physical, performance, and technical data

Physical characteristics	
Weight	260 g (9.2 oz) including batteries
Dimensions	Length: 15.7 cm Width: 5.5 cm Height: 3.6 cm
Technical performance data	
Max. pressure (pump)	150 cmH2O
Pressure set range	5 cmH2O to 50 cmH2O
Default pressure	25 cmH2O
Resolution (setting/display)	±1 cmH2O
Pressure accuracy	±2 cmH2O
Electrical specifications	
AC power input	100 to 240 VAC, 50 to 60 Hz ± 10% 1.25 VA typical, 3.25 VA maximum
USB input power	5 VDC ± 5% 1.0 W typical, 2.5 W maximum
Car adapter	12 V to 24 V ± 5% 1.25 W typical, 3.25 W maximum
Fuse	A fuse is integrated in the adapter. Replace only with same type and value (T 3.15 A, L 250V, 5x20).
Battery requirements	AA (IEC-HR6) NiMH rechargeable 1.2 V, ≥1900 mA
Alarm volume ¹	55.0 dB(A) ± 6 dB(A)
1. Volume measured at 1 m distance from the IntelliCuff device	
Environmental conditions	
Relative humidity	Storage: 5% to 95%, noncondensing Operation: 5% to 95%, noncondensing

Temperature	<p>Operation: -15°C to 50°C (5°F to 122°F)</p> <p>NOTE. If the the device was stored at temperatures exceeding 50°C (122°F), cool it down for at least 15 minutes in an ambient temperature of 20°C (68°F) before use.</p> <p>Storage: -15°C to 70°C (5°F to 158°F)</p>
Altitude	<p>-650 to 7620 m (-2132 to 25,000 ft) above sea level</p> <p>Above 4000 m (13,123 ft) only battery power is supported.</p>
Noise	< 30 dB(A)

	Indicates the degree of protection against electric shock according to IEC 60601-1. Class II devices have double or reinforced insulation, as they have no provision for protective grounding
	The TÜV NRTL mark with the indicators “C” and “US” means that the product complies with Canadian requirements and the requirements of US authorities for safety
	Dispose according to Council Directive 2002/96/EC or WEEE (Waste Electrical and Electronic Equipment)
	Serial number
	Reference number

12.2 Symbols on the device labels

Symbol	Description
	Manufacturer
	Date of manufacture
	IntelliCuff poses unacceptable risks to patient, medical staff, or other persons within the MR environment.
	Consult operator’s manual. Refer to the operator’s manual for complete information. This label on the device points the user to the operator’s manual for complete information. In the operator’s manual, this symbol cross-references the label.
	Symbol for “Caution”
	Applicable to pediatric/adult patient groups
	Applicable to pediatric patient group
CE 0197	CE Marking of Conformity, seal of approval guaranteeing that the device is in conformance with the Council Directive 93/42/EEC concerning medical devices
	Type BF applied part (classification of medical electrical equipment, type BF, as specified by IEC 60601-1)

	Device is limited to operating temperatures of -15°C to 50°C (5°F to 122°F)
	Recyclable materials
IP 34	Protected against splashing water and solid particles larger than 2.5 mm
	Do not use any blades, knives, or cutters to open; they can damage the product

Symbols on package for PN 282016

Symbol	Description
	Not made with natural rubber latex
	Do not reuse
	Do not use if the packaging is damaged
	Batch code
	Quantity
	Use-by date

12.3 System overview

The IntelliCuff device is an automatic cuff pressure controller. It provides the operator with continuous monitoring and the means to adjust cuffed tracheal tubes and cuffed tracheostomy tubes.

12.4 Essential performance

The applied cuff pressure must be maintained and monitored. If it is higher or lower than the set limits, this must be detected and the operator informed through an alarm.

13 Parts and accessories

Name	Part number
Car adapter	PN 951011
Cuff pressure tube (2 m) with filter	PN 282016
Hanger for IntelliCuff	PN 951020
Rail mount	PN 951055
Multi mount	PN 951057
Click mount adapter	PN 951019
Battery cover	PN 951016
Power supply kit	PN 951012

14 Warranty

LIMITED WARRANTY

THE WARRANTY DESCRIBED IN THIS AGREEMENT IS IN LIEU OF ANY AND ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. HOWEVER, IMPLIED WARRANTIES ARE NOT DISCLAIMED DURING THE PERIOD OF THIS LIMITED WARRANTY.

Hamilton Medical guarantees its products to be shipped free from defects in material and workmanship. The warranty does not include disposable items. Disposable items and consumable products are considered to be of single use or of limited use only and must be replaced regularly as required for proper operation of the product following the operator's manual.

Hamilton Medical and the manufacturer shall have no obligations nor liabilities in connection with the product other than what is specified herein, including without limitation, obligations and/ or liabilities for alleged negligence, or for strict

liability. In no event shall the company be liable for incidental or consequential damages, either direct or contingent.

This Limited Warranty shall be void and not apply:

1. If the product has not been installed and connected by an authorized local representative of Hamilton Medical in accordance with the instructions furnished by Hamilton Medical and by a Hamilton Medical representative.
2. If no evidence is present that the occurrence of damage/ repair happened within the certified warranty period.
3. If the serial number has been altered, effaced or removed and there is no bill of sale or evidence to verify the product's purchase date.
4. If the defects arise from misuse, negligence, or accidents or from repair, adjustment, modification or replacement made outside Hamilton Medical's factories or other than an authorized service center or authorized service representative.
5. If the product has been modified, or in any nature altered without prior written authorization from Hamilton Medical.
6. If the product is or has been used in any way that is not specified under "Intended Use".
7. If the product has been used by anyone but properly trained personnel under the supervision of a physician.

Replacements and/or repairs furnished under this Limited Warranty do not carry a new warranty, but carry only the unexpired portion of the original Limited Warranty. The warranty of repaired and/or replaced components does not exceed the Limited Warranty of the device.

To obtain service under this Limited Warranty, claimant must promptly notify the country's sales partner of Hamilton Medical regarding the nature of the problem, serial number and the date of purchase of the Product.

Except as stated above, Hamilton Medical shall not be liable for any damages, claims or liabilities including, but not limited to, personal bodily injury, or incidental, consequential, or special damages. Nor will Hamilton Medical be liable for any damages, claims or liabilities including, but not limited to, personal bodily injury, or incidental, consequential, or special damages resulting from misuse of the device or failure to comply with any of the provisions made in this manual.

14.1 Miscellaneous

The general terms and conditions of Hamilton Medical shall be applicable. This agreement shall be governed by and construed in accordance with the laws of Switzerland and may be enforced by either party under the jurisdiction of the court of Chur, Switzerland.

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Hamilton Medical, Inc.