

On the move. HAMILTON-HF90 for patient mobilization

Two hangers for water bottle and infusion bags

Sturdy handles to support your patient and maneuver the trolley

Adjustable height for clip-in trolley mount



Easy-to-maneuver wheels with stoppers



Adjustable telescope pole with clip to attach and store the circuit

Two standard rails for additional accessories or support arm

Cylinder holder for up to three oxygen bottles



On the move. With our optional and interchangeable battery

Oxygen consumption on the main screen

Battery status on the main screen

Up to two hours of battery runtime¹

¹ With the following settings: Adult/Ped patient group, Flow = 40 l/min, Oxygen = 30%, Temperature = 37°C, SpO2 sensor connected.

Battery status indicator on the battery



For a personal demonstration, contact your account manager or email us at info@hamilton-medical.com!

EL020250103N.00

The products shown here are not available for purchase by the general public. The information provided here is intended for healthcare professionals only. Always read the labels and follow the product's instructions for use. Specifications are subject to change without notice. Some features are options. Not all features are available in all markets. All images are used for illustrative purposes only and may not accurately represent the product or its use. For all proprietary trademarks (®) and third-party trademarks used by Hamilton Medical AG see www.hamilton-medical.com/trademarks. ©2025 Hamilton Medical AG. All rights reserved.